


I-Series Straight Double cuff


VERIFIED



Series	Internal segment	Cuffs	Length	Medcomp Kit	A (cm)	B (cm)	C (cm)	D (cm)	E (cm)	J (cm)	K (cm)
I-series	Straight (Tenckoff)	Double	31 cm	MDC31S	6.5	0.8	3	0.8	20	3	31
I-series	Straight (Tenckoff)	Double	32.5 cm	MDC32S	7	0.8	3.5	0.8	20	5	32.5
I-series	Straight (Tenckoff)	Double	37 cm	MDC37S	10	0.8	5	0.8	20	5	37
I-series	Straight (Tenckoff)	Double	42 cm	MDC42S MPD-242	15	0.8	5	0.8	20	7	42
I-series	Straight (Tenckoff)	Double	47 cm	MDC47S	20	0.8	5	0.8	20	15	47

I-Series Coiled Double cuff

VERIFIED



Series	Internal segment	Cuffs	Length	Medcomp Kit	A (cm)	B (cm)	C (cm)	D (cm)	E (cm)	F (cm)	G (cm)	K (cm)	I (cm)
I-series	Coiled (Tenckoff)	Double	57.5 cm	MDC57C MPD-257	12	0.8	6	0.8	18	5	8	57.5	31
I-series	Coiled (Tenckoff)	Double	63 cm	MDC63C MPD-263	18	0.8	6	0.8	17	5	14	63	38

I-Series Straight Single cuff (subcutaneous)




Series	Internal segment	Cuffs	Length	Medcomp Kit	A (cm)	D (cm)	E (cm)	J (cm)	K (cm)
I-series	Straight (Tenckoff)	Single (Subcutaneous)	31 cm	MSC31S	10	0.8	20	3	31
I-series	Straight (Tenckoff)	Single (Subcutaneous)	37 cm	MSC37S	16	0.8	20	5	37
I-series	Straight (Tenckoff)	Single (Subcutaneous)	42 cm	MSC42S	21	0.8	20	7	42
I-series	Straight (Tenckoff)	Single (Subcutaneous)	46 cm	MSC46S MPD-146S	25	0.8	20	15	46

VERIFIED


I-Series Straight Single cuff (Preperitoneal)

VERIFIED



Series	Internal segment	Cuffs	Length	Medcomp Kit	A (cm)	D (cm)	E (cm)	J (cm)	K (cm)
I-series	Straight (Tenckoff)	Single (Preperitoneal)	41 cm	MPP41S MPD-141	15	0.8	25	7	41
I-series	Straight (Tenckoff)	Single (Preperitoneal)	46 cm	MPP46S MPD-146	20	0.8	25	15	46

I-Series Coiled Single cuff (Pre-peritoneal)



Series	Internal segment	Cuffs	Length	Medcomp Kit	A (cm)	D (cm)	E (cm)	F (cm)	G (cm)	I (cm)	K (cm)
I-series	Coiled (Tenckoff)	Single (Preperitoneal)	40.25 cm	MPP40C	16	0.8	12	4	13	26	40
I-series	Coiled (Tenckoff)	Single (Preperitoneal)	57.5 cm	MPP57C	12	0.8	26	5	7	32	58
I-series	Coiled (Tenckoff)	Single (Preperitoneal)	60 cm	MPP60C MPD-160	15	0.8	25	5	10	34	60