




Silicone	Coiled	15F
Material	Tip Design	French Size



COILED TIP

LARGE SIDE HOLES

TITANIUM BARBED CONNECTOR



CATHETER ONLY		
MC20XC	15F X 64CM X-SERIES, COILED	5/BOX

BASIC SET		
MC20XS	15F X 64CM X-SERIES, COILED	5/BOX

ACCESSORIES		
MP217	17F VASCU-SHEATH® INTRODUCER SET	5/BOX
MC50	MALE LUER LOCK END CAP	100/BOX
MBCA	BARBED LUER LOCK CATHETER ADAPTOR	20/BOX
MCS	CATHETER STYLET	5/BOX
MTS	TUNNELING STYLET	10/BOX
MTC	TITANIUM LUER ADAPTOR	5/BOX
TBC	TITANIUM BARBED CONNECTOR	50/BOX
TRX	X-SERIES TUNNELING ROD	5/BOX
FSX	X-SERIES FALLER STYLET	5/BOX
40501	PREOPERATIVE MARKING STENCIL PACK	1/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Titanium Barbed Connector
- (1) Female Luer
- (1) Clamp
- (1) End Cap

BASIC SET CONTAINS:

- (1) Catheter
- (1) Titanium Barbed Connector
- (1) Female Luer
- (1) Clamp
- (1) End Cap
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Tunneler
- (1) 10cc Syringe
- (6) 4" x 4" Gauze
- (1) Tape Measure

Refer to the Table of Contents for Important Risk Information regarding this device.